

Bano toilet grab bars / back rest / shower rail / wall grab bars / shower seats / accessories

Care & maintenance instructions



BANO AS

Industrivegen 22
6823 Sandane-Norway
Phone: +47 57 86 98 00
post@bano.no

Bano toilet grab bars - Back rest

Care & maintenance instructions

Products:

5600/R-02 right and **5600/L-02** left 900 mm

Bano grey toilet grab bar- height adjustable

5607/R-02 right and **5607/L-02** left 700 mm

Bano grey toilet grab bar- height adjustable

5600/R-B-02 right and **5600/L-B-02** left 900 mm

Bano black grab bar- height adjustable

5607/R-B-02 right and **5607/L-B-02** left 700 mm

Bano black toilet grab bar- height adjustable

5601-02 900 mm Bano grey toilet grab bar- non adjustable

5608-02 700 mm Bano grey toilet grab bar- non adjustable

5601-B-02 900 mm Bano black toilet grab bar- non adjustable

5608-B-02 700 mm Bano black toilet grab bar- non adjustable

5603 Bano toilet grab bar- 900 mm- grey

For h.adj toilet/cistern 5970-EL and 5970-MA

5604 Bano toilet grab bar- 700 mm- grey

For h.adj toilet/cistern 5970-EL and 5970-MA

5603-B Bano toilet grab bar- 900 mm- black

For h.adj toilet/cistern 5970-EL and 5970-MA

5604-B Bano toilet grab bar- 700 mm- black

For h.adj toilet/cistern 5970-EL and 5970-MA

5603-EL Bano toilet grab bar with eletrical height control - 900 mm- grey

For toilet/cistern 5970-EL and 5970-MA

5604-EL Bano toilet grab bar with eletrical height control - 700 mm- grey

For toilet/cistern 5970-EL and 5970-MA

5972-04 Bano back rest - grey

5972-04-B Bano back rest - black



5600/L-02



5601-02



5603



5972-04



Material specification

- **Toilet grab bars:** structural foam with polyurethane coating
- **Fixing rails:** - natural anodised aluminium alloy
- galvanized steel
- **Plastic caps:** injection-moulded polypropylene
- **Back rest:** structural foam with polyurethane coating

Maintenance instructions

Standard cleaning agents can be used for cleaning. Use a soft cloth or sponge with a mild cleaning agent. Spray onto the unit and wipe off with a moist cloth. Do not use scouring powder or steel wool, as this will damage the finish. Cleaning agents that contain alcohol and are used in the pharmaceutical industry can be used.

Warning

Never use scouring agents or cleaning agents that contain strong chemicals. Never use acids, methanol, methylated spirits, white spirit or lynol.

Thank you for taking the time to read these instructions carefully – it is time well spent!

Products:

5581/L-G Bano shower rail/grab bar- left- grey
5581/R-G Bano shower rail/grab bar- right- grey

5581/L Bano shower rail/grab bar- left- white
5581/R Bano shower rail/grab bar- right- white

5581/L-B Bano shower rail/grab bar- left- black
5581/R-B Bano shower rail/grab bar- right- black


5591/ L-G Bano wall grab bar L shape- left- grey
5591/ R-G Bano wall grab bar L shape- right- grey
5591/ L Bano wall grab bar L shape- left- white
5591/ R Bano wall grab bar L shape- right- white
5591/ L-B Bano wall grab bar L shape- left- black
5591/ R-B Bano wall grab bar L shape- right- black

5440-G Bano wall grab bar 400 mm- grey
5460-G Bano wall grab bar 600 mm- grey
5480-G Bano wall grab bar 800 mm- grey
54100-G Bano wall grab bar 1000 mm- grey

5440 Bano wall grab bar 400 mm- white
5460 Bano wall grab bar 600 mm- white
5480 Bano wall grab bar 800 mm- white
54100 Bano wall grab bar 1000 mm- white

5440-B Bano wall grab bar 400 mm- black
5460-B Bano wall grab bar 600 mm- black
5480-B Bano wall grab bar 800 mm- black
54100-B Bano wall grab bar 1000 mm- black

5311-G Bano grab bar for cabinets- grey
5311-W Bano grab bar for cabinets- white
5311-B Bano grab bar for cabinets- black



Material specification

• **Shower rail and wall grab bars:**

Anti-bacterial, powder-coated aluminium, nylon and stainless steel screws

• **Caps:** nylon, ABS plastic

Maintenance instructions

Standard cleaning agents can be used for cleaning. Use a soft cloth or sponge with a mild cleaning agent. Spray onto the object and wipe off with a moist cloth. Do not use scouring powder or steel wool, as this will damage the finish. Cleaning agents that contain alcohol and are used in the pharmaceutical industry can be used.

Warning

Never use scouring agents or cleaning agents that contain strong chemicals. Never use acids, methanol, methylated spirits, white spirit or lynol.

Thank you for taking the time to read these instructions carefully – it is time well spent!

Products:

- 5700-02** Bano shower seat with height adjustment- grey
- 5700-B-02** Bano shower seat with height adjustment- black
- 5701-02** Bano shower seat without height adjustment- grey
- 5701-B-02** Bano shower seat without height adjustment- black
- 5710-02** Bano armrest- grey
- 5710-B-02** Bano armrest- black



Material specification

- **Seat and armrest:** structural foam with polyurethane coating
- **Fixing rails:** natural anodised aluminium alloy
- **Plastic caps:** injection-moulded polypropylene

Maintenance instructions

Standard cleaning agents can be used for cleaning. Use a soft cloth or sponge with a mild cleaning agent. Spray onto the unit and wipe off with a moist cloth. Do not use scouring powder or steel wool, as this will damage the finish. Cleaning agents that contain alcohol and are used in the pharmaceutical industry can be used.

Warning

Never use scouring agents or cleaning agents that contain strong chemicals. Never use acids, methanol, methylated spirits, white spirit or lynol.

Thank you for taking the time to read these instructions carefully – it is time well spent!

Bano accessories

Care & maintenance instructions

Products:

3477-G Bano wall mounted white toilet brush- grey holder

3477 Bano wall mounted white toilet brush- white holder

5400-G Bano single wall hook- grey

5400 Bano single wall hook- white

5400-B Bano single wall hook- black

5402-G Bano double wall hooks- grey

5402 Bano double wall hooks- white

5402-B Bano double wall hooks- black

5491-G Bano bathroom shelf- grey

5491 Bano bathroom shelf- white

5491-B Bano bathroom shelf- black

5596-G Bano shower shelf- grey

5596 Bano shower shelf- white

5596-B Bano shower shelf- black

5496-G Bano waste bin 10L- grey

5496 Bano waste bin 10L- white

5510 Bano walking stick holder

3477-G



5402-G




5400-G



5596-G




5491-G



5496-G



5510



Material specification

- **Shower self / bathroom self:** Nylon
- **Single/ double hooks:** Cast in solid nylon
- **Waste bin:** Nylon
- **Toilet brush:** Polypropylene plastic **Holder:** ABS-plastic
- **Walking stick holder:** ABS plastic and rubber

Maintenance instructions

Standard cleaning agents can be used for cleaning. Use a soft cloth or sponge with a mild cleaning agent. Spray onto the object and wipe off with a moist cloth. Do not use scouring powder or steel wool, as this will damage the finish. Cleaning agents that contain alcohol and are used in the health care service can be used.

Warning

Never use scouring agents or cleaning agents that contain strong chemicals. Never use acids, methanol, methylated spirits, white spirit or lynol.

Thank you for taking the time to read these instructions carefully – it is a time well spent!